

With the same features and benefits as the OmniHolder Family, the OmniHolder-MINI is an indispensable accessory to improve order, patient safety, and infection control.



Why use The **OmniHolder-Mini**

The OmniHolder-Mini eliminates potential medical complications from mispositioning the pressure transducers leading to significant errors in the measurement of invasive pressures that will result in incorrect therapeutic decisions and increased morbidity and mortality. By securing the OmniHolder-Mini at the correct site of the ICU bedrail, the pressure transducers are always calibrated at heart level, regardless of the bed position.

The OmniHolder-Mini contributes to infection control by eliminating the traditional "custom-made roll towels" that hold the transducers with tape and which often get soiled by bodily fluids. By eliminating the "roll towels," there are fewer objects and materials on the bed touching the patient. It opens bed space for more clarity in patient care in critical moments.

The OmniHolder-Mini provides ergonomics and order by having all the pressure lines secure in grommets and the infusion drugs organized on a lateral manifold for better visualization of the IV lines. This feature contributes to patient safety by eliminating potential errors of injecting fluids or drugs in the wrong line.

Who can **benefit**



Hospitals









Patients

Physicians Nurses

Insurance



Top Hook
Self-adjusting rotational
mechanism adapting to any
back incline up to 45 degrees.



Brackets
Two clip dockets to secure an arterial and venous transducer.



GrommetsTo secure and organize pressure lines




Manifolds pockets
Two detachable lateral pockets
to secure multiple drips.





The OmniHolder-Mini



Technical specifications:

Dimensions	Weight capacity
Height: 163 mm (6.4") Width: 121 mm (4.7")	8.5 lb

Ordering info and support. We will help you